

## B 6.- Miscellaneous Accessories

Miscellaneous Accessories consist of a variety of aftermarket handles to fit on most any power wheelchair's joystick. Each selection enables users of different abilities access to power mobility.

### B 6.1.- BioForm Joystick Handles

#### B 6.1.1.- BioForm Goal Post Joystick Handles

BioForm Goal Post joystick handles ( 'Y' shape), provides a superb option for greater control of a power wheelchair for the individual who lacks dexterity with their hands. They are shaped from a piece of 1/16" x 1" aluminum with a coupling piece on the underside for easy attaching to the joystick. This coupling device has a set screw for securing it to the joystick.




**BioForm Goal Post Joystick Handle**

Model	Description	Size	Dimension		
			W	GH	H
JGS	Goal Post Joystick Handle w/set screw	Small	3	1	2.5
JGM	Goal Post Joystick Handle w/set screw	Medium	4	1	3
JGL	Goal Post Joystick Handle w/set screw	Large	5	1	3.5
JGX	Custom Goal Post Joystick Handle	Custom	Template or Specs		

#### B 6.1.2.- BioForm 'T' Joystick Handles

BioForm 'T' joystick handles provide an option for excellent control of a power wheelchair for the individual who lacks dexterity with their hands. They are shaped from a 1/2" aluminum rod with a coupling piece on the underside for easy attaching to the joystick. This coupling device has a set screw for securing it to the joystick.



**BioForm 'T' Joystick Handle**

Model	Description	Size	Dimension		
			W	D	H
JTS	'T' Joystick Handle w/set screw, 3" wide	Small	3	0.5	2.5
JTM	'T' Joystick Handle w/set screw, 4" wide	Medium	4	0.5	3
JTL	'T' Joystick Handle w/set screw, 5" wide	Large	5	0.5	3.5
JTX	Custom 'T' Joystick Handle	Custom	Template or Specs		

***"If you don't see it.... Contact Us"***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.